

SUOMEN LAITOSTEKSTIILIT KY

Asemakatu 25

FI-94100 KEMI, FINLAND

Tel. +358 16 282 580

Fax +358 16 282 580

GSM: +358 40 7778000

E-mail: suomen.laitostekstiilit@co.inet.fi

Mattress covers:

Material: Polyurethane-coated fabric. Non-sweating, waterproof, feels pleasant, does not swish, wash at 95 °C, can be wiped, hard-wearing. **FIREPROOF, ANTIMILDEW, ANTIFUNGAL.**

Also very suitable as a replacement cover of a hygiene mattress.

All sizes available.

Type PS1

With an elastic band



Type PS4

With a zipper



Pillowcase:

Type TS1

Material: See mattress covers

With a flap



Rescue sheet:

Type EL-2A

Placed under the mattress.

Quick and easy to take into use,
and light to trail.

Small elastic bands keep the Velcro tapes
neatly in place.

Material: Slippery nylon with PU coating.

Wash at 80 °C, do not squeeze dry.

Size: 195 x 75 cm



Transfer mattress:

Type SP1

Used for transferring a patient from a bed
to another bed or to a shower stretcher.

Easy to use thanks to its slipperiness.

In addition, the material is waterproof so that
it can easily be wiped clean.

Material: Slippery nylon with polyurethane
coating + 2-cm-thick expanded foam inside.

Wash at 80 °C, do not squeeze dry.

Size: 180 x 60 cm



Instructions for use: Two persons on opposite sides of the bed. One turns the patient to his/her side and at the same time lifts one side of the bed cover, turn sheet etc. up. Simultaneously, the other pushes the transfer mattress as far as possible under the patient so that the patient's head is on the transfer mattress. Lower the patient down. By holding to the bed cover, one person pulls the patient onto the carriage next to the bed. The other one can lift the other end of the bed cover. The patient lightly slides onto the carriage. Finally, remove the transfer mattress.

Apron:

Type ES1

Material: See mattress covers

Non-sweating, very suitable for e.g. bath attendants.



Food covers:

Type RS1

Material: See mattress covers

Size: 45 x 90 cm

Fastened with poppers

The edge of the bag has a plastic strip to
keep the bag continuously open.

Type RS2

Same as RS1, but the bag is in the middle.



Multi-purpose belt:

With padding that can be wiped clean
Soft, comfortable in use. Quick-release fastener.
Material: See mattress covers + 2-cm-thick expanded foam



Separate leg section

Prevents the patient from sliding off the chair.
Material: See mattress covers + 0.5-cm-thick expanded foam

Sling:

Type LSK1 for hands

Soft and non-sweating, easy to wipe clean.
Material: See mattress covers + 0.5-cm-thick expanded foam
Width 10 cm



Type LSJ1 for legs

Width 15 cm

Lifting sheet for shower use:

Type NS1

Material: Nylon with PU coating.
Easy to keep clean.
Wash at 80 °C, do not squeeze dry.
Size: 100 x 100 cm



Lifting mat:

Type NL1

Handy aid for lifting a patient out of a chair.

